**lorelli**<sup>®</sup>  
designed in EU

*Non-Contact IR Thermometer*

*For Body & Surface*

V1.6W

## MANUAL INSTRUCTION

find us on  
[www.lorelli.eu](http://www.lorelli.eu)

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# IMPORTANT! KEEP FOR FUTURE REFERENCE! READ CAREFULLY!

EN

## INTRODUCTION

**Thank you for buying this non-contact forehead infrared thermometer.** Please read very carefully this instruction before using the device because the correct use is of critical significance for the accurate measurement of temperature. This thermometer is carefully designed for accurate, safe, and fast measurement of the body temperature on the forehead. Following the guidance in this instruction will make the thermometer user-friendly. The forehead temperature is the same as the subcutaneous blood flow. This is good for checking the body temperature because measuring is not compromised by external factors and delays like measuring by the oral and rectal method.



## PRECAUTIONS AND SAFETY MEASURES

1. Sweat on the forehead may result in incorrect values. Before measurement, please ensure that the forehead is clean. If you have cleaned the forehead, wait 5 - 10 minutes before measuring the temperature.
2. Do not touch or blow over the infrared sensor. The contamination of the infrared sensor may result in incorrect readings. Clean the sensor if you think that it might be contaminated.  
(see section "Care and cleaning")
3. To achieve correct measurements, it is important to keep the probe clean and without scratches. Finger stains and dirt will affect the accuracy of the thermometer. To achieve accurate measurements, keep the probe's tip clean. After cleaning, place the thermometer in a dry place at room temperature. Do not expose the thermometer to extreme temperatures, moisture, direct sunlight. Protect it from dropping and impacts.
4. **ATTENTION!** Do not service and repair the device while in use.
5. This device is not life-supporting.
6. The use of this thermometer does not replace the professional consultation with your family doctor. Measurements used only for reference. Consult your doctor if you have any doubts about your health status.
7. If the ambient temperature in the storage room considerably differs from the ambient temperature at the place where measurement will be made, wait a while so that the temperature of the thermometer becomes equal to the ambient temperature. In case of big differences in the ambient temperatures, condition the thermometer for about 30 minutes before use.
8. Appropriate to be used at home conditions.
9. Keep the thermometer away from children and pets to avoid ingestion or inhalation of small parts. Do not allow children to measure temperature unattended. Children may not be capable to use the device in compliance with the instruction for use. This thermometer is not a toy.
10. Do not store the thermometer under extreme temperatures (below -4°F /-20°C and above 131°F/55°C) or under extreme moisture (Rh>95% moisture condensation). Please also observe the operating conditions: the following temperatures 10°C - 40°C (50°F-104°F) Humidity: ≤85%RH, moisture condensation; Atmospheric pressure: 70-106 kPa. The thermometer under these extreme conditions may be inaccurate.
11. Do not use the thermometer if it is damaged in any way (for example - a broken infrared sensor). Continuous use of a damaged device may cause injury, incorrect readings, or serious hazards.
12. Do not disassemble, repair or replace any parts of the device except change of batteries.
13. Do not use a mobile phone near the thermometer while making measurements.
14. If the thermometer will not be used for a longer period, remove the batteries to prevent possible damage caused by chemical leakage. In case of leakage from the battery, carefully remove them.
15. In case of any leakage from the batteries, protect your skin. If this happens during use, do not touch the patient and the batteries simultaneously.

16. This device will not cause allergic reactions or irritation during normal use. However, if you experience any allergic reactions to any material, stop using the thermometer and consult your doctor.
17. Store the thermometer in a dry and clean place. Do not expose the product to any chemical solvents, dust, moss, direct sunlight, or high temperatures.
18. Do not store the thermometer together with objects with sharp edges.
19. Do not touch wounds, skin affected by any skin disease or trauma with the sensor of the thermometer.
20. It is not advisable that the thermometer is used by patients that have been treated with anti-inflammatory drugs for skin, to place the sensor of the thermometer on skin exposed to direct sunlight, heated by a source of heat, exposed to direct flow of an air-conditioner and by patients treated with a cold compress.
21. If the following conditions are available, it is recommended that the temperature is measured at least three times and to use the highest measured value:
- newborn babies under three months.
  - children under the age of three years with an impaired immune system.
  - if you use the infrared thermometer for the first time.
22. The thermometer is not appropriate for sterilization. Do not use it in OXYGEN SATURATED ENVIRONMENT!
23. Please notify us of any unexpected situations.
24. After dropping, impact, or other circumstances that may affect the readings, please refer to a service centre.
25. This device shall be serviced, repaired, and opened only at authorized service centres.
26. This thermometer is not intended for preterm babies.

## PRODUCT DESCRIPTION

**Principle of operation of the thermometer** - the temperature sensor collects the infrared energy emitted by the skin surface. Once it is focused by lens, energy is transformed into a temperature reading by the measuring chains.

**Body temperature** - the normal body temperature varies in certain limits. The trend is that the normal temperature decreases with ageing. The table below shows the normal temperature ranges depending on age.

0 – 2 years	36.4 – 38.0 °C	97.5 – 100.4 °F
3 – 10 years	36.1 – 37.8 °C	97.0 – 100.0 °F
11 – 65 years	35.9 – 37.6 °C	96.6 – 99.7 °F
> 65 years	35.8 – 37.5 °C	96.4 – 99.5 °F

The normal temperature range varies among people and can be influenced by many factors such as the time of the day, measuring at different places, level of activity, drugs, emotions, etc.


We advise you to practice measuring the temperature on yourself and the family members when in good health. Thus you will learn how to use the thermometer and will feel more confident when measuring the temperature of the same people when they feel sick.

**Features** - Non-contact design, safe, and more hygienic to use. Fast measuring, for less than 1 second. Accurate and reliable. User-friendly, single-button design. Multifunctional, capable of measuring forehead, room, milk, water, object temperature. Memory function for 35 values, easy to retrieve. Switching between unmute and mute mode. Function for high temperature, shown in orange and red light. Switching between °C and °F. Automatic switch off and energy saving.

**Measuring the forehead temperature** - press the button for measuring to switch on the thermometer directed to the forehead at a distance of 0 - 3 cm. Skin contact is not necessary. When the button for measuring the body temperature is pressed, a sound is issued and then you may see the reading on the display. Holding the thermometer too long in the hand before measuring may result in the warming of the device. It means that measuring can be incorrect.

**NOTE:** The forehead measurement is indicative. The forehead temperature may vary upwards to 1 ° F / 0,5 ° C from your actual body temperature. Please consider the factors that affect accuracy as it is described in the instruction. If the eyebrows area is covered with hair, sweat, or dirt, please clean it to improve the measuring accuracy. Always ensure that the sensor is clean. Always ensure that the user and the thermometer were in the same room before measuring.


### Product description


1. LCD
2. Measuring button
3. Memory button
4. Mute button (or °C/°F switching)
5. Sensor
6. Batteries cover

### Description of the display

1. Temperature value
2. Object temperature mode
3. Forehead temperature mode
4. Fahrenheit/ Celsius degrees
5. Battery level
6. Recalling values from memory
7. Unmute/ mute

**Measuring object temperature** - When the thermometer is switched off, press the Mem button for 3 seconds. Then press the room or object measuring button. Keep the thermometer at a distance of 1-3 cm from the object. Press and release the measuring button for 1 second and when you hear the sound signal, you may see the reading.

**Readings** - „T“ shows measuring temperature in forehead mode.

1. If  $32^{\circ}\text{C} \leq T \leq 37.3^{\circ}\text{C}$  ( $89.6^{\circ}\text{F} \leq T \leq 99.2^{\circ}\text{F}$ ), the green light will stay for 12 seconds, with a single long sound signal.
2. If  $37.4^{\circ}\text{C} \leq T \leq 37.9^{\circ}\text{C}$  ( $99.3^{\circ}\text{F} \leq T \leq 100.3^{\circ}\text{F}$ ), the orange light will stay for 12 seconds, with 3 short sound signals and the value on the LCD will blink, which is a warning that you have a mild fever.
3. If  $38^{\circ}\text{C} \leq T \leq 42.9^{\circ}\text{C}$  ( $100.4^{\circ}\text{F} \leq T \leq 109.2^{\circ}\text{F}$ ), the red light will stay for 12 seconds, with 5 short sound signals and the value on the LCD will blink, which is a warning that you have a high temperature.

**Switching between unmute and mute mode** - when the thermometer is switched on, press the mute button and then again to unmute.

**Check of the 35 data sets from the memory** - when the thermometer is switched on, press the Mem button to enter the memory mode, press again that button to check the 35 data sets one by one. If no values are available, it will show "--- M".

**°C / °F conversion** - when the thermometer is switched on, the continuous pressing of the mute button for 5 seconds will switch °C / °F.

**Temperature compensation setting** - when the thermometer is switched on, press simultaneously the mute button and the Mem button for 2-3 seconds to enter temperature compensation mode. Temperature adjustment from  $\pm 0,0$  to  $\pm 2,0$  is made by pressing the Mem button.

Note: All future temperatures, which you will measure, will be automatically added to the value that you have set.

**Switching off** - the device will be automatically switched off when idle for 10 seconds. Or by pressing the measuring button continuously for 6 seconds.

**Attention:**

1. All values from the memory will be reset when you remove the batteries.
2. All settings will be set to default when you remove batteries. If you have to adjust the settings, please switch on and make the new settings.

**Changing batteries** - slide the cover of the batteries in the indicated direction. Place correctly two batteries AAA in the department. Remove the batteries if you will not use the thermometer for more than two months.

**Attention:** Under certain conditions, batteries may leak and this may cause chemical burning or damage the device. If the leakage from the batteries contacts the skin or the eyes, immediately flush with plenty of water. Batteries must be handled under the supervision of an adult. Dispose of the used batteries with due care. Do not disassemble, recharge or throw batteries in a fire. The battery may explode with flames. The use of alkaline batteries is recommended. Keep the batteries out of the reach of children. Do not mix old and new batteries or batteries of a different type.

**Disposal and recycling** - Batteries contain substances that are harmful to the environment. Do not throw the batteries together with the ordinary household waste but deliver them to a specialized place for collecting batteries. Always remove the battery before disposal of the device or delivering it to an official waste collection point. **IMPORTANT!** After the end of the service life of the device, do not dispose of it together with the ordinary household waste but deliver it to an official collection point for recycling. Thus you will help with environmental protection.

### **Care and cleaning**

1. Cleaning and user maintenance must not be performed by children.
2. Remove the batteries from the device before cleaning.
3. Use a pad soaked with spirit or cotton pad, soaked in 95% alcohol to clean the case of the thermometer and the measuring sensor.
4. Once the alcohol fully dries, you can measure again (after min. 10 minutes).
5. Ensure that no liquid enters inside the thermometer. If this happens, use the device again once it has fully dried out.
6. Do not use abrasive cleaning agents, solvents, or benzene for cleaning, and never soak the device in water or other cleaning liquids.
7. Be careful and avoid scratching the surface of the LCD.



**Guarantee and after-sales service** - all damages caused by incorrect use as a result of non-observance of the instructions in the operations manual will not be covered by the guarantee.

### **Electromagnetic compatibility (EMC)**

1. This product needs special precautions concerning EMC and must be installed and operated in compliance with the provided information about EMC and this device can be affected by portable and mobile radio communication equipment.
2. Do not use mobile telephones or other devices, which create electromagnetic fields near the device. This may cause faulty operation of the device.
3. Attention: This device has been thoroughly tested and inspected to ensure its correct functioning and operation!
4. Attention: This thermometer shall not be used next to or in contact with other equipment. If the equipment must remain in immediate proximity to the used thermometer, it shall be monitored to confirm the normal operation during measurement.

By increasing the number of electronic devices, such as computers and mobile telephones, the used medical devices may be affected by electronic disturbances from other devices. Electromagnetic disturbances may affect the correct operation of the medical device and create a potentially hazardous situation. Medical devices also shall not interfere with other devices. To regulate the EMC (electromagnetic compatibility) requirements to avoid situations with hazardous products, EN 60602-1-2 standard has been introduced. This standard determines the levels of protection against electromagnetic disturbances, as well as the maximum levels of electromagnetic emissions for medical devices. This medical device has been produced and complies with the standard EN60601-1-2 on protection and emissions.

Symptom	Possible Cause	Description & Solution
Failed to power on.	The battery level is too low.	Replace with a new battery
	Polarities of the batteries are reversed.	Ensure the batteries are in the right position
	The thermometer is damaged	Contact dealer
The reading is too low	The lens of the probe is dirty.	Clean the lens with a cotton swab.
	The distance of the item and target is too far	Ensure the thermometer measuring the forehead center within 3cm.
	You have just come from a cold environment	Stay in a warmer room for at least 30 minutes before taking a reading
The reading is too high	You have just come from a hot environment.	Stay in an adequately cool room for at least 30 minutes before taking a reading
	The ambient temperature is not in range.	3 short beeps and red backlit for 3 seconds. Take a measurement under an ambient temperature between 10°C (50°F) and 40°C (104°F).
	Memory Error	3 short beeps and red backlit for 3 seconds. Contact dealer.
	In forehead mode, $T > 42.9^{\circ}\text{C}$ (109.2°F)	3 short beeps and red backlit for 3 seconds.
	In forehead mode, $T < 32^{\circ}\text{C}$ (89.6°F)	3 short beeps and red backlit for 3 seconds.
	$2.5\text{V} \pm 3\% \leq \text{power voltage} \leq 2.6\text{V} \pm 3\%$	The battery level is low, it suggests you to replace the battery, but you can continue to use it.
	The power voltage is lower than $2.5\text{V} \pm 3\%$ .	It will turn off automatically after 30 seconds. Please replace with a new battery

## Specifications

Product name	Infrared thermometer	
Power supply	DC1.5V×2	
Measurement range	Forehead: 32.0°C–42.9°C (89.6°F–109.2°F)	
	Object: 0°C–100°C (32°F–212°F)	
Accuracy (Laboratory)	Forehead mode	±0.2°C /±0.4°F
	Object mode	±1.0°C/1.8°F
Display resolution	0.1°C/°F	
Measuring distance	0-3cm /0-1.5 Inches	
Automatic shutdown	10s±1s	
Memory	35 groups of measured temperature.	
Operational conditions	Operating environment: 10°C - 40°C (50°F-104°F) ≤85%RH Relative Humidity: ≤85%RH, moisture condensation Atmospheric pressure: 70-106 kPa	
Storage and transport conditions	RH ≤95% -20°C ~ 55°C(-14°F~131°F)Atmospheric pressure:70-106 kPa	
Battery	2*AAA, can be used for more than 3000 times	
Weight & Dimension	69 g (without battery),157x35x36 mm	
Device Life	5 Years	

## Symbol

Symbol	Description
	Type BF applied part.
	Information about a manufacturer
	Please read the instructions carefully.
IP22	2 Protected against solid foreign objects of 12,5 mm Ø and greater
	<b>IMPORTANT</b> Inaccurate reading or thermometer damage may occur if the thermometer is not correctly used.
	Waste electrical materials should be sent to a dedicated collection point for recycling.
	EU safety requirements
	Radio frequency emitter symbol
	European authorized representative
SN	Serial number
LOT	Batch number

## **IEC 60601-1-2 ME EQUIPMENT and ME SYSTEMS identification, marking and documents for Class B product**

EMC information

IEC 60601-1-2 ME EQUIPMENT and ME SYSTEMS identification, marking and documents for Class B product

The ME EQUIPMENT or ME SYSTEM is suitable for home healthcare environments and so on.

**Warning:** Don't near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.

**Warning:** Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

**Warning:** Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.”

**Warning:** Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Infra-red Thermometer (FC-IR200), including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

**If any:** a list of all cables and maximum lengths of cables (if applicable), transducers and other ACCESSORIES that are replaceable by the RESPONSIBLE ORGANIZATION and that are likely to affect compliance of the ME EQUIPMENT or ME SYSTEM with the requirements of Clause 7 (EMISSIONS) and Clause 8 (IMMUNITY). ACCESSORIES may be specified either generically (e.g. shielded cable, load impedance) or specifically (e.g. by MANUFACTURER and EQUIPMENT OR TYPE REFERENCE).

**If any:** the performance of the ME EQUIPMENT or ME SYSTEM that was determined to be ESSENTIAL PERFORMANCE and a description of what the OPERATOR can expect if the ESSENTIAL PERFORMANCE is lost or degraded due to EM DISTURBANCES (the defined term “ESSENTIAL PERFORMANCE” need not be used).

## Technical description

1.all necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to electromagnetic disturbances for the excepted service life.

2. Guidance and manufacturer's declaration -electromagnetic emissions and Immunity

**Table 1**

Guidance and manufacturer's declaration - electromagnetic emissions	
Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	N/A
Voltage fluctuations/ flicker emissions IEC 61000-3-3	N/A

**Table 2**

Guidance and manufacturer's declaration - electromagnetic Immunity		
Immunity Test	IEC 60601-1-2 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Electrical fast transient/burst IEC 61000-4-4	Power supply lines: 2kV input/output lines: ±1 kV 100 kHz repetition frequency	N/A
Surge IEC 61000-4-5	line(s) to line(s): ±1 kV. line(s) to earth: ±2 kV.	N/A
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% 1 cycle And 70% 25/30 cycles Single phase: at 0 0% 250 cycle	N/A
Power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC61000-4-6	150KHz to 80MHz: 3Vrms 6Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	N/A
Radiated RF IEC61000-4-3	10 V/m 80 MHz – 2,7 Ghz 80 % AM at 1 kHz	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz
NOTE U <sub>i</sub> is the a.c. mains voltage prior to application of the test level.		

**Table 3**

Guidance and manufacturer's declaration - electromagnetic Immunity							
	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Modulation (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	385	380 –390	TETRA 400	Pulse modulation 18Hz	1,8	0.3	27
	450	430–470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28
	710	704 – 787	LTE Band 13, 17	Pulse modulation 217 Hz	0,2	0.3	9
	745						
	780						
	810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
	870						
	930						
	1720	1 700–1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
	1845						
	1970						
	2450	2 400 –2 570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
	5240	5 100 –5 800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0,2	0.3	9
	5500						
	5785						

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Product Service

## EU Quality Management System Certificate (MDR)

Pursuant to Regulation (EU) 2017/745 on Medical Devices, Annex IX Chapters I and III  
(Class IIa and Class IIb Devices)

**No. G10 004665 0011 Rev. 00**

**Manufacturer:** **Shenzhen Finicare Co., Ltd**  
201, No. 50, the 3rd Industrial Park  
Houting Community, Shajing Street  
Bao'an District  
518104 Shenzhen  
PEOPLE'S REPUBLIC OF CHINA

SRN Manufacturer - CN-MF-000024364

**Authorized  
Representative:**

Riomavix S. L.  
Calle de Almansa 55, 1D, 28039 Madrid, SPAIN

The Certification Body of TÜV SÜD Product Service GmbH certifies that the manufacturer has established, documented and implemented a quality management system as described in Article 10 (9) of the Regulation (EU) 2017/745 on medical devices. Details on device categories covered by the quality management system are described on the following page(s). The Report referenced below summarises the result of the assessment and includes reference to relevant CS, harmonized standards and test reports. The conformity assessment has been carried out according to Annex IX Chapter I and III of this regulation with a positive result. The quality management system assessment was accompanied by the assessment of technical documentation for devices selected on a representative basis. The certified quality management system is subject to periodical surveillance by TÜV SÜD Product Service GmbH. The surveillance assessment shall also include an assessment of the technical documentation for the device or devices concerned on the basis of further representative samples. All applicable requirements of the testing and certification regulation of TÜV SÜD Group have to be complied with.

For details and certificate validity see: [www.tuvsud.com/ps-cert?q=cert:G10\\_004665\\_0011\\_Rev\\_00](http://www.tuvsud.com/ps-cert?q=cert:G10_004665_0011_Rev_00)

**Report No.:** GZ2336501

**Valid from:** 2023-10-17

**Valid until:** 2028-10-16

Christoph Dicks  
Head of Certification/Notified Body

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Фирма : Shenzhen Finicare Co., Ltd, 201, No.50, the 3rd Industrial Park, Houting Community, Shajing Street, Bao'an District, Shenzhen 518104, China

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For and on behalf of  
SHENZHEN FINICARE CO., LTD.

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## Declaration of Compliance

Company Name : Shenzhen Finicare Co., Ltd, 201, No.50, the 3rd Industrial Park, Houting Community, Shajing Street, Bao'an District, Shenzhen 518104, China

Product Item No: FC-IR200

Description : Infrared Thermometer

For our model FC-IR200, delivered to DIDIS Ltd as 1025014 NON-CONTACT IR THERMOMETER FOR BODY AND SURFACE:

1. We manufacturer hereby declare that the above model is compliant to the Restriction of Hazardous Substances in Electrical and Electronic Equipment Directive (RoHS 2.0) 2011/65/EU.
2. We hereby agree that only RoHS compliant materials and components will be used in the assembly of the above model. Under no circumstances will non-compliant materials and components be used in place of approved compliant equivalents during the life of the product.
3. This device corresponds to the Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 on medical devices, amending Directive 2001/83/EC, Regulation (EC) No 178/2002 and Regulation (EC) No 1223/2009 and repealing Council Directives 90/385/EEC and 93/42/EEC relating to medical devices, ISO 80601-2-56 and European standard EN 60601-1-2 and it is subject to specific precautions, regarding the electromagnetic compatibility.

Signature and stamp: .....

*For and on behalf of*  
SHENZHEN FINICARE CO., LTD.

Position: ..... *Chau Li* .....

Date: 19<sup>th</sup>, December 2023 ..... *Authorized Signature(s)*

Manufacturer / Производител:  
Shenzhen Finicare Co., Ltd  
201, No.50, the 3rd Industrial Park,  
Houting Community, Shajing Street,  
Bao'an District, Shenzhen 518104  
China

Share Info Consultant  
RIOMAVIX SO0CIEDAD LIMITADA  
Calle de Almansa 55, 1D, Madrid, 2839 Spain  
ES-AR-000001202

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#### **Didis Ltd.**