# Tender Care

Thermometer

# Digital

## MANUAL INSTRUCTION

designed in EU

[www.lorelli.eu](http://www.lorelli.eu)

V1.3W

find us on

<b>EN</b>	Manual Instruction.....	3
<b>BG</b>	Инструкция за употреба.....	7
<b>HU</b>	Használati utasítás.....	11
<b>NL</b>	Gebruikshandleiding.....	15
<b>ES</b>	Instrucciones de uso.....	19
<b>IT</b>	Istruzione per l'uso.....	23
<b>FR</b>	Mode d'emploi.....	27
<b>GR</b>	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ ΚΑΡΕΚΛΑΣ.....	31
<b>AR</b>	دليل لتيلا تعليمات.....	35
<b>DE</b>	Bedienungsanleitung.....	43
<b>CZ</b>	Návod k použití.....	47
<b>RU</b>	Инструкция по эксплуатации.....	51
<b>SRB HR ME BIH</b>	Uputstvo za upotrebu.....	55
<b>RO</b>	Instrucțiuni de utilizare.....	59
<b>AL</b>	Insruksion për përdorim.....	63
<b>TR</b>	Kullanım talimatı.....	67
<b>PL</b>	Instrukcja użytkownika.....	71
<b>MK</b>	Упатствата за употреба.....	75

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**OPERATING INSTRUCTIONS**

PLEASE READ THESE INSTRUCTIONS CAREFULLY BEFORE USING

The digital thermometer offers a safe, accurate and quick temperature reading. You can measure your temperature through oral or underarm. Thoroughly read this instruction manual before using the digital thermometer. This thermometer is manufactured according to the system of EN 80601-2-56 and complies with the requirements of EN60601-1, EN 60601-1-2. The manufacturer of this digital thermometer is ISO13485 and MDD 93/42/EEC certified

**INFORMATION OF PREDICTIVE OFFSET**

Although the generally accepted "normal" temperature is 37.0°C (98.6°F), temperatures can vary from 36.1°C (96.9°F) to 37.2°C (98.9°F) and still may be considered "normal". Variations in temperature can be attributed to activities such as exercise, smoking, eating and drinking. Even time of day may influence your temperature. For example, your temperature is lower in the morning than in the afternoon. Other variations may be due to the location of the temperature reading. While the oral temperature follows the guidelines described above, an axillary (under the arm) temperature will be 0.5°C (1.0°F) lower.

**PRODUCT DESCRIPTION****IMPORTANT SAFETY NOTES**

To assure the correct use of the product basic safety, measures should always be followed including the precautions listed below.

**WARNING**

- A high or prolonged fever requires medical attention, especially for young children. Please contact your doctor.
- Carefully read and follow the enclosed instructions to ensure accurate temperature readings. Note that temperature readings are affected by many factors including physical exertion, the drinking of hot or cold beverages before measurement, as well as measuring technique.
- Please keep still during measurement.
- The use of temperature readings for self-diagnosis is dangerous. Consult with your doctor for the interpretation of results. Self diagnosis may lead to the worsening of existing disease conditions.
- This thermometer is used for taking temperatures through oral or underarm. Do not attempt to take temperatures at other sites, such as in the ear, as it may result in false readings and may lead to injury.
- Store the thermometer out of the reach of children. Do not allow children to take their temperatures unattended. Children may injure themselves when attempting to take temperatures without supervision.
- Do not leave the battery, battery cap or probe cover where children can get to them. Children may swallow them. Should a child swallow the battery, battery cap or probe cover, contact a doctor immediately.
- Do not attempt measurements when the thermometer is wet as inaccurate readings may result.
- Do not repair or maintain the thermometer when in use
- Please put a special trash can to recycle used batteries, and send them to a facility that specializes in disposing of used batteries. Do not discard them at will, otherwise they will pollute the environment and water sources.
- Measurement results are affected by damage to probe or insufficient battery power
- Do not place in high temperature, direct sunlight, and contact with any chemical solvent to prevent chemical changes and affect operation function
- It takes at least 30 minutes for the product to get to normal working conditions at high or low temperatures

**CAUTION**

- Do not bite on the thermometer. Doing so may lead to breakage and/or injury
- Do not share the thermometer among individuals.
- Do not attempt to disassemble or repair the thermometer. Doing so may result in inaccurate readings.
- Please do not repair the thermometer by yourself. Please contact the manufacturer or the manufacturer's representative when the thermometer reports an error and requires maintenance or needs guidance on setting operations.
- Do not attempt to incinerate the battery. It may burst.
- Pay attention to polarity (+ -) when replacing the battery. Failure to do so may lead to fluid leakage, heat generation or bursting, resulting in damage to the unit.
- Remove the battery when the thermometer will not be used for 3 months or more. Failure to do so may lead to fluid leakage, heat generation or bursting, resulting in damage to the unit.
- Do not use mobile phones near the thermometer.
- Do not use the thermometer in places where strong static electricity or electromagnetic fields are present. Doing so may lead to inaccurate readings and may contribute to instrument failure.
- Do not step on the unit or the protective case.
- Do not attempt to disinfect the sensing section of the thermometer by immersing in alcohol or in hot water (water over 50°C).
- Contraindication: It is prohibited for infants and young children, mental disorders, coma, oral disorders, nose and mouth surgery, mouth breathing; trauma, surgery, inflammation under the armpit, a lot of sweating under the armpit, shoulder joint injury or weight loss, the thermometer is not tightly clamped to use the thermometer.
- Degradation of sensor performance may result in inaccurate measurement. Please contact the manufacturer.

## GENERAL SAFETY PRECAUTIONS

- The patient is the intended operator, Do not use the thermometer to measure anything other than human body temperatures.
- Do not drop the thermometer and protective case or expose them to shock or vibration.
- Do not store the thermometer in the protective case when it is damp. Wipe it off with a dry cloth first.
- Avoid taking the temperature until after 30 minutes has lapsed after exercise, bathing or eating/drinking.

## SWITCHING BETWEEN FAHRENHEIT/CELSIUS

In the off state, press and hold the power button for 3 seconds to switch between Fahrenheit and Celsius

### HOW TO USE

1. Press ON/OFF button to activate. The unit will beep and display. This LCD display test will go on for about 2 seconds.
2. When Lo and flashing °C(°F) displays, the thermometer is ready for temperature measuring.
3. If the room temperature is higher than 32.0°C(89.6°F), the room temperature will be displayed instead of Lo°C(Lo°F).
4. The beep will sound when temperature measuring is complete. The degree sign of °C(°F) on the LCD will stop flashing.
5. During measurement: LCD will display "Lo°C" or "Lo°F". If the measured temperature is below 32.0°C(89.6°F); The unit will display "Hi°C" or "Hi°F" if the measured temperature is above 42.9°C/109.2°F.
6. The thermometer will automatically turn off in about 10 minutes for energy-saving. It is suggested to press the On/Off button to turn the thermometer off once the temperature measuring is finished.

## SPECIFICATIONS

Model	T15SC
Range:	32.0°C-42.9°C (89.6°F-109.2°F)
Accuracy:	±0.1°C, 35.5°C-42.0°C (±0.2°F, 95.9°F-107.6°F) ±0.2°C under 35.5°C or over 42.0°C (±0.4°F under 95.9°F or over 107.6°F) at standard room temperature of 25°C (77.0°F)
Display:	Liquid crystal display
measurement time	≤1min
Battery life:	About 1000 times (22°C)
Memory:	For storing the last measured value
Battery:	DC 1.5 V (size LR41)
Power consumption:	0.15 mW (in measurement mode)
Dimension:	12.4cm x 1.8cm x 0.9cm(L x W x H)
Weight:	Approx. 10 grams including battery
Using environment:	Temperature: 10°C-40°C (50°F-104°F) Relative humidity: 15%RH-85%RH Atmospheric pressure: 86kPa-106kPa
Storage and transportation condition:	Temperature: -25°C -55°C (-13°F-131°F) Relative humidity: 15%RH-85%RH Atmospheric pressure: 50kPa-106kPa
Type of protection against electric shock	Internal powered equipment.
Degree of protection against electric shock	Type BF applied part.
Classification according to the degree of protection against ingress of water	Ip22( The first number 2: Protected against solid foreign objects of 12,5 mm Φ and greater. The second number: Protected against vertically falling water drops when enclosure tilted up to 15°.)
Equipment not suitable for use in the presence of flammable mixtures	Not AP or APG equipment
Mode of operation:	Continuous
Software version	T11_V2.0
Applied part	Digital Thermometer surface
Intended apply human body parts	Oral and AXILLARY

## ORAL USE

The mouth should remain closed up to 2 minutes before attempting a reading.

\* Place the probe tip in the mouth under the tongue so that it rests to the left or right of the root of the tongue.

\* Use downward tongue pressure to hold the thermometer in place.

\* Hold the thermometer to keep it from sliding around in the mouth.

The normal body temperature measured by this method is  $36.3^{\circ}\text{C}$ - $37.2^{\circ}\text{C}$

( $97.3^{\circ}\text{F}$ - $98.9^{\circ}\text{F}$ ) Approximate measuring time: 1minute.

## AXILLARY USE

Wipe the underarm with a dry towel and close your armpit at least 5 minutes before attempting a reading.

\* Place the probe tip under the arm so the tip is touching the skin with the thermometer perpendicular to the body.

Position the arm across the chest so the probe tip is well covered by the arm.

\* The normal body temperature measured by this method is  $36.0^{\circ}\text{C}$ - $37.0^{\circ}\text{C}$  ( $96.8^{\circ}\text{F}$ - $98.6^{\circ}\text{F}$ )

Approximate measuring time: 1minute.

Note: Axillary temperatures are generally  $0.5^{\circ}\text{C}$ / $1.0^{\circ}\text{F}$  lower than oral readings.

## EXPLANATION OF SAFETY SIGNS AND SYMBOLS

\* Type BF applied part

\* Symbol for "THE OPERATION GUIDE MUST BE READ"

Symbol for "MANUFACTURER"

Symbol for "COMPLIES WITH MDD 93/42/ECCREQUIREMENTS"

Cautions: Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the equipment or other property


Warnings: Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury

DISPOSAL: Do not dispose this product as unsorted municipal waste. Collection of such waste separately for special treatment is necessary. \* OPERATOR do not touch the battery and the patient simultaneously

## BATTERY REPLACEMENT

1. When the  appears in the lower right corner of the LCD, the battery is exhausted and need to be replaced.
2. Pull the battery cap off with hand in the direction shown below.
3. Gently pull out the circuit board with battery compartment approximately 1 cm.
4. Remove the depleted battery with a pointed object such as a pen. Discard used batteries out of reach of children.
5. Place a new 1.5V DC button size battery type LR41 or equivalent in the compartment with positive pole side faced up and negative pole side faced down
6. Close the battery cover.

### Attention!

Lid, sealants, batteries are small parts, children should be noted that upon contact, do not inhale or swallow.

### Warning!

No modification of this equipment is allowed

## CLEANING INSTRUCTIONS

Before and after each use, clean the thermometer with a soft cloth and isopropanol diluted with water. Do not immerse the thermometer in fluids or sterilize with boiling method, gas or steam autoclave.

It is necessary to clean and disinfect the same unit between different users.

## MAINTENANCE

If you need help, please contact us or our representative for help.

## Handling of common problems

Use exception, if first of all, please check the following points

phenomenon	reason	resolvent
Turn on the power Screen no display	Battery exhaustion	Replacement of new batteries
	Battery polarity is not correct	Battery polarity consistent with battery holder identification
Low measurement temperature	Wrong measurement position	Measure according to instructions
	Temperature measuring probe or dirt in temperature measuring area	Please measure after removing dirt
ERR	Damaged temperature probe	Return to manufacturer for maintenance
Lo Hi	When measured $\leq 31.9^{\circ}\text{C}$	Please read the instruction manual re-measure them
	When the measurement $\geq 43.0^{\circ}\text{C}$	

GUANGDONG GENIAL TECHNOLOGY CO., LIMITED  
**GUANGDONG GENIAL  
 TECHNOLOGY CO., LIMITED**

1-6-05-02, 11th Road, Area B,  
 Guangfozhao Economic  
 Cooperation Zone,  
 Zhagang Town, Huaiji County,  
 Zhaoqing City, Guangdong  
 Province, 526437, P.R. China

Table 1

## CONTENT OF EMC

1\* WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally."

2\* WARNING: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation."

3\* WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the digital thermometer, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

declaration - electromagnetic emission	
Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B

Table 2

declaration - electromagnetic immunity		
Immunity test	IEC 60601 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	$\pm 8$ kV contact $\pm 2$ kV, $\pm 4$ kV, $\pm 8$ kV, $\pm 15$ kV air	$\pm 8$ kV contact $\pm 2$ kV, $\pm 4$ kV, $\pm 8$ kV, $\pm 15$ kV air
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m

NOTE: UT is the a.c. mains voltage prior to application of the test level.

Table 3

declaration - electromagnetic immunity		
Immunity test	IEC 60601 test level	Compliance level
Radiated RF IEC 61000-4-3	10V/m 80 MHz to 2.7 GHz	10V/m

Table 4

declaration - IMMUNITY to proximity fields from RF wireless communications equipment					
Immunity test	IEC60601 test level				Compliance level
	Test frequency	Modulation	Maximum power	Immunity level	
Radiated RF IEC 61000-4-3	385 MHz	**Pulse Modulation: 18Hz	1.8W	27 V/m	27 V/m
	450 MHz	*FM+ 5Hz deviation: 1kHz sine	2 W	28 V/m	28 V/m
	710 MHz	**Pulse Modulation:	0.2 W	9 V/m	9 V/m
	745 MHz	217Hz			
	780 MHz	217Hz			
	810 MHz	**Pulse Modulation:	2 W	28 V/m	28 V/m
	870 MHz	18Hz			
	930 MHz	18Hz			
	1720 MHz	**Pulse Modulation:	2 W	28 V/m	28 V/m
	1845 MHz	217Hz			
1970 MHz	217Hz				
2450 MHz	**Pulse Modulation: 217Hz	2 W	28 V/m	28 V/m	28 V/m
5240 MHz	**Pulse Modulation:	0.2 W	9 V/m	9 V/m	9 V/m
5500 MHz	217Hz				
5785 MHz	217Hz				

Note\* - As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.  
 Note\*\* - The carrier shall be modulated using a 50 % duty cycle square wave signal.

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MedPath GmbH
**Didis Ltd.**